



## 🮾 SOL-M™ Blunt Fill Needle





## **SOL**-M<sup>™</sup> Blunt Fill Needle with 5 micron filter



| Ref                                      | Size | Gauge | Length             | Qty/Box | Qty/Case |
|--|------|-------|--------------------|---------|----------|
| 110022F                                  | —    | 18    | 1 <sup>1</sup> /2" | 100     | 1000     |
| Blunt Fill Needle with Luer Lock Syringe |      |       |                    |         |          |
| 1831815B                                 | 3mL  | 18    | 1 1/2"             | 100     | 800      |
| 1851815B                                 | 5mL  | 18    | 1 1/2"             | 100     | 800      |
| 1811815B                                 | 10mL | 18    | 1 <sup>1</sup> /2" | 100     | 1200     |

Sol-M<sup>™</sup> Blunt Fill Needles are used for draw-up of medication, re-constitution of medication and accessing IV Bags.



Sol-M<sup>™</sup> Blunt Fill Needle with filter features a 5 micron filter mem-brane that is designed to prevent larger particles being drawn into the syringe when preparing medications from ampules.

## Features & Benefits

- 45° bevel to help reduce accidental needlesticks
- 18G for ease in draw-up of medications
- 5 micron filter Effective filtration of glass particles
- Available individually or pre-attached to 3ml, 5ml and 10ml syringes

